

We make our patients' orthotics with Riecken's PQ Gel.

ORTHOTIC
LABORATORY

Make custom orthotics with Riecken's Professional Quality Viscoelastic Polymer Gel

PQ Liquid is available in 2 part, pint or quart kits and mixes easily for consistent results.

- A true, soft-tissue supplement
- Ideal for arthritics and diabetics
- Make total contact, custom inlays and/or flexible orthotics, distal end pads and bases for healing sandals
- Make custom pads any durometer
- Cures in approximately 20 minutes at 72°
- Contains no nutritive components for microbe growth
- Contains no cellular constructions so there is no collapse in load-bearing areas
- Fast and simple to use when custom casting molds
- Requires no expensive production equipment.

PQ produces a custom inlay with maximum shock absorption and dissipation. Its uses include flexible orthotics, distal end pads and healing sandals.

View the Ecatalog of

- Insoles
- Heel Cups
- Pads
- Lifts
- Wedges
- Spring Plates

at www.FootComfortWorld.net

Carmela Anders
C. Ped.



Call (800) 331-8040, ext 100, to order PQ Gel Foot-Care Products.

Riecken's PQ
Laboratory of Foot-Care Products
Evansville, IN

PQ Gel Foot-Care Products are also available from Cascade, G&W Heel Lifts, PEL, P.W. Minor, R&B Medical, Ruby Leather, Verne Bintz and Warwick Enterprises.